

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human C1QTNF6 Polyclonal Antibody

多克隆抗体

产品货号: PA12834

第三版

描述: This protein is not uniformly distributed but rather focally localized in hepatocellular carcinoma cells and the tumor endothelial cells. It may play a role in metabolism and immunity. CTRPs acts as new molecular mediators connecting inflammatory and metabolic diseases.

抗原: Recombinant protein of human C1QTNF6

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: