

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human VWF Polyclonal Antibody

多克隆抗体

产品货号: PA12853

第三版

描述:The glycoprotein encoded by this gene functions as both an antihemophilic factor carrier and a platelet-vessel wall mediator in the blood coagulation system. It is crucial to the hemostasis process. Mutations in this gene or deficiencies in this protein result in von Willebrand's disease. An unprocessed pseudogene has been found on chromosome 22.

抗原:Recombinant protein of human VWF

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Plasma.