

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human A4GNT Polyclonal Antibody

多克隆抗体

产品货号: PA12862

第三版

描述:This gene encodes a protein from the glycosyltransferase 32 family. The enzyme catalyzes the transfer of N-acetylglucosamine (GlcNAc) to core 2 branched O-glycans. It forms a unique glycan, GlcNAc α 1 \rightarrow 4Gal β 1 \rightarrow 6R and is largely associated with the Golgi apparatus membrane.

抗原:Recombinant protein of human A4GNT

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in stomach and pancreas.