
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ARHGEF5 Polyclonal Antibody

多克隆抗体

产品货号: PA12920

第三版

描述:Rho GTPases play a fundamental role in numerous cellular processes initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein may be involved in the control of cytoskeletal organization.

抗原:Recombinant protein of human ARHGEF5

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Mainly expressed in kidney, liver, pancreas, lung and placenta.