

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ANKRD54 Polyclonal Antibody

多克隆抗体

产品货号: PA12956

第三版

**描述:** Plays an important role in regulating intracellular signaling events associated with erythroid terminal differentiation. Interacts (via ankyrin repeat region) with LYN (via SH3-domain) in an activation-independent status of LYN. By similarity. Forms a multiprotein complex with LYN and HCLS1. By similarity. Interacts with TSN2, VAV1, DBNL AND LASP1.

**抗原:** Recombinant protein of human ANKRD54

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**