

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

4-1BB ligand, 71-254aa, Human, His-tagged, Recombinant, E.coli

产品货号: TP01299

第三版

别名:Tumor necrosis factor ligand superfamily member 9(TNFSF9); 4-1BBL; CD137L.

描述:4-1BB ligand is a type 2 transmembrane glycoprotein belonging to the TNF superfamily. 4-1BBL is expressed by activated B cells, macrophages, dendritic cells, activated T cells, neurons and astrocytes.

配方:Liquid. In 20 mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 20% glycerol

分子量:23.8 kDa (222aa)

序列:

MRGSHHHHHHGMASMTGGQQMGRDLYDDDDKDRWGSHMREGPELSPDDPAGLLDLRQGMFAQLVAQNVLLIDGP
LSWYSDPGLAGVSLTGGLSYKEDTKELVVAKAGVYYVFFQLELRRVVAGEGSGSVSLALHLQPLRSAAGAAALALTVDLPP
ASSEARNSAFGFQGRLLHLSAGQRLGVHLHTEARARHAWQLTQGATVLGLFRVTPEIPAGLPSRSE

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.