



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

EMCN, 19-190aa, Human, His tag, E.coli,

产品货号: TP02013

第三版

别名:Endomucin, EMCN2, MUC14.

描述:Endomucin, also called endothelial sialomucin or mucin-like sialoglycoprotein, which interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix. EMCN is a mucin-like sialoglycoprotein that interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix.

配方:Liquid. In 20 mM Tris-HCl buffer, pH8.0, 10% glycerol, 1mM DTT, 100mM NaCl

分子量:20.5kDa (197aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMNSTGVLEAANNSLVVTTTKPSITTPNTESLQKNVVPTTGTTPKGTITNELLKM
SLMSTATFLTSKDEGLKATTTDVRKNDSIISNVTVTSTLPLNAVSTLQSSKPKTETQSSIKTTEIPGSVLQPDASPSKTGTLTSIP
VTIPENTSQSQVIGTEGGKNASTSATRSYSS

纯度:> 95% by HPLC

浓度:0.25mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.