



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

FABP5, 1-135aa, Human, Recombinant, E.coli

产品货号: TP02076

第三版

别名:Fatty acid binding protein 5, EFABP, E-FABP, PAFABP, PA-FABP

描述:FABP5 is a member of the intracellular fatty acid binding protein (FABP) family, which is known for the ability to specifically bind fatty acids (FAs) with high affinity for stearic and linoleic acids. FABP5 is expressed in endothelial cells of the microvasculature of the placenta, heart, skeletal muscle, small intestine, lung, and renal medulla.

配方:Liquid in 20 mM Tris pH 8.0, 1 mM DTT, 20 % Glycerol

分子量:15 kDa (135 aa)

序列:

纯度:> 95% by HPLC

浓度:1 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.