

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

FAIM, 1-213aa, Human, His tag, E.coli

产品货号: TP02083

第三版

**别名:** fas apoptotic inhibitory molecule 1 isoform a, FAIM1

**描述:** FAIM protects against death receptor-triggered apoptosis and regulates B-cell signaling and differentiation. Several transcript variants encoding different isoforms have been found for this gene. Recombinant human FAIM protein, fused to His-tag at N-terminus, was expressed in E.coli.

**配方:** Liquid. 20mM Tris-HCl buffer (pH8.0) containing 10% glycerol 0.4M Urea

**分子量:** 26.4kDa (236aa)

**序列:**

MGSSHHHHHHSSGLVPRGSHMGSM LLPFIRTLPLLCYNHLLVSPDSATLSPPY SLEKMTDLVAVVDVALSDGVHKIEFEH  
GTTSGKR VVYVDGKEEIRKEWMFKLVGKETFYVGA AKTKATINIDAI SGFAYEY TLEINGKSLKKY MEDRSKTTNTWVLH  
MDGENFRIVLEKDAMDVWCNGKKLETAGEFVDDGTETHFSIGNHDCYIKAVSSGKRKEGIIHTLIVDNREIPEIAS

**纯度:** > 95% by HPLC

**浓度:** 0.5mg/ml (determined by Bradford assay)

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.