



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

FAM107B, 1-131aa, Human, His tag, E.coli

产品货号: TP02084

第三版

别名:FAM107B protein, FLJ45505, MGC11034

描述:FAM107B protein, also known as FAM107B, is expressed in most tissues in the human body without there being a high frequency in any one tissue. It is found in all stages of human development. Recombinant human FAM107B protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

分子量:17.9 kDa (154aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSM AEPDYIEDDNPELIRPQKLINPVKTSR NHQDLHRELLMNQKRGLAPQNKPELQK
VMEKRKRQDQVIKQKEEEA QKKKSDLEIELLKRQQKLEQLELEKQKLQEEQENAPEFVKVKGNLRRTGQEVAAQES

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.