

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Gankyrin, 1-226aa, Human, Recombinant, E.coli

产品货号: TP02194

---

第三版

别名: proteasome 26S subunit non ATPase 10 (PSMD10) isoform 1, ankyrin repeat protein, p28

描述: Gankyrin is an oncoprotein commonly overexpressed in most hepatocellular carcinomas. Gankyrin interacts with S6 ATPase of the 19S regulatory particle of the 26S proteasome. It is known to play a role in the regulation of the phosphorylation of the retinoblastoma protein by CDK4, and to enhance the ubiquitinylation of p53 by MDM2.

配方: Liquid in Phosphate buffered saline (pH7.4), 10% glycerol

分子量: 24.4kDa (226aa)

序列:

纯度: > 95% by HPLC

浓度: 1 mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.