

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

GTSF1, 1-167aa, Human, His tag, E.coli

产品货号: TP02374

第三版

别名: Gametocyte-specific factor 1, FAM112B

描述: GTSF1 belongs to the UPF0224 (FAM112) family and contains 1 CHHC-type zinc finger. Recombinant human GTSF1 protein fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT.

分子量: 21.7 kDa (190aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMEEETYTDSLDPKLLQCPYDKNHQIRACRFPYHLIKCRKNHPDVASKLATCPFNAR
HQVPRAEISHHISSCDDRSCIEQDVVNQTRSLRQETLAESTWQCPCDEDWDKDLWEQTSTPFVWGTTHYSDNNSPASN
IVTEHKNNLASGMRVPKSLPYVLPWKNNNGNAQ

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.