

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

HAVCR2, 22-202aa, Human, His tag, E.coli

产品货号: TP02405

第三版

别名:Hepatitis A virus cellular receptor 2, FLJ14428, KIM-3, Tim-3, TIM3, TIMD3

描述:HAVCR2, also known as TIM3, belongs to the immunoglobulin superfamily. HACVR2 regulates macrophage activation. It Inhibits T-helper type 1 lymphocyte (Th1)-mediated auto- and alloimmune responses and promotes immunological tolerance. Recombinant human HAVCR2 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

分子量:22.7 kDa (206aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMSEVEYRAEVGQNAYLPCFYTPAAPGNLVPVCWGKACPVFECGNVVLRTDER
DVNYWTSRYWLNDFRKGDVSLTIENVTLADSGIYCCRIQIPGIMNDEKFNLLKVIKPAKVTPAPTLQRDFTAAFPRMLTT
RGHGPAETQTLGSLPDINLTQISTLANELRDSRLANDLRDSGATIRIG

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.