



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

HDGF, 1-100aa, Human, E.coli

产品货号: TP02425

第三版

别名:Hepatoma-derived growth factor isoform a, HMG1L2

描述:Hepatoma-derived growth factor isoform a, also known as HDGF, belongs to the hepatoma-derived growth factor family. HDGF was originally identified as a secreted mitogen from the human hepatoma cell line Huh- 7. This nuclear targeted vascular smooth muscle cell mitogen is a heparin-binding protein that is highly expressed in tumor cells where it stimulates proliferation.

配方:Liquid. In 20mM Tris-HCl buffer (pH7.5) containing 1mM DTT,10% glycerol, 1mM EDTA

分子量:11.5 kDa (100aa)

序列:

MSRSNRQKEYKCGDLVFAKMKGYPHWPARIDEMPEAAVKSTANKYQVFFFGTHETAFLGPKDLFPYEEESKEKFGKPNKR  
KGFSEGLWEIENNPTVKASGY

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.