

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

HDGFL1, 1-251aa Human, His tag, E.coli

产品货号: TP02426

第三版

别名:Hepatoma-derived growth factor-like protein 1, DJ309H15.1, PWWP1

描述:HDGFL1 belongs to the HDGF family and contains 1 PWWP domain. Recombinant human HDGFL1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid, In Phosphate buffered saline (pH7.4) containing 10% glycerol

分子量:29.6 kDa (274aa) confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMSAYGMPMYKSGDLVFAKLKGYAHWPARIEHMTQPNRYQVFFFGTHETAFLSPK  
RLFPYKECKEKFQKPNKRRGFSAGLWEIENNPTVQASDCPLASEKSGDGPWPEPEAAEGDEDKPTHAGGGGDELGKPD  
DDKPTEEEKGPLKRSAGDPPEDAPKRPKEAAPDQEEEEAEAEAAEAERAAAAAATAVDEESPFLVAVENGSAPSEPLVC  
EPPQPEEEELREEEVADEEASQEWHAEPGGGDRDSL

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.