

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Hexokinase 3 Human, His-tagged, Recombinant, E.coli

产品货号: TP02443

第三版

别名: Hexokinase type III, HK III

描述: Hexokinase is the first enzyme in the glycolytic pathway, catalyzing the transfer of a phosphoryl group from ATP to glucose to form glucose-6-phosphate and ADP. In mammals, four distinct enzymes-types 1 to 4 hexokinases-have been identified. The enzyme is found in most cells, but there is tissue specificity for the particular type of hexokinase. Hexokinase3 lacks the hydrophobic N-terminal sequence critical for targeting to mitochondria.

配方: Liquid. In 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

分子量: 101.1kDa (943aa)

序列:

```
MGSSHHHHHHSSGLVPRGSHMDSIGSSGLRQGEETLSCSEELPGPSDSELVQECLQQFKVTRAQLQQIQASLLGSMEQA  
LRGQASPAPAVRMLPTYVVGSTPHGTEQGDFVLELGATGASLRVLWVTLTGIEGHRVEPRSQEFVIPQEVMLGAGQQLF  
DFAAHCLSEFLDAQPVNKQGLQLGFSFSPCHQTGLDRSTLISWTKGFRCSGVEGQDVVQLLRDAIRRQGAYNIDVVAV  
VNDTVGTMMSGCEPGRVPCVGLVVDGTNACYMEEARHVAVLDEDRGRVCSVEWGSFSDDGALGPVLTTFDHTLD  
HESLNPGAQRFEKMIGGLYLGLVRLVLAHLARCVLFGGCTSPALLSQGSILLEHVAEMEDPSTGAARVHAILQDLGLSP  
GASDVELVQHVC AAVCTRAAQLCAAALAAVLSCLQHSREQQTLQVAVATGGRVCERHPRFCSVLQGTVMMLAPECDV  
SLIPVDGGGRGVAMVTAVAARLAAHRRLEETLAPFRLNHDQLAAVQAQMRKAMAKGLRGEASSLRMLPTFVRATPD  
GSERGDFLALDLGGTNFRVLLVRVTTGVQITSEIYSIPETVAQGSQQQLFDHIVDCIVDFQKQGLSGQSLPLGFTFSFPCR  
QLGLDQGILLNWTGFKASDCEGQDVVSLREAITRRQAVELNVVAIVNDTVGTMMSGYEDPRCEIGLIVGTGTNACY  
MEELRN VAGVPGDSGRMCINMEWGAFGDDGSLAMLSTRFDASVDQASINPGKQRFKEMISGMYLGEIVRHILLHLTSLG  
VLFRRGQQIQRLQTRDIFKTKFLSEIESDSLALRQVRAILEDLGLPLTSSDALMVLEVCQAVSQRAAQLCGAGVAAVVEKIRE  
NRGLEELAVSVGVDGTYKLVHPRFSSLVAATVRELAPRCVVTFLOQSEDGSGKGAALVTAVACRLAQLTRV
```

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.