

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

HRAS-like suppressor3,1-133aa, Human,Recombinant, E.coli

产品货号: TP02501

第三版

别名:HRASLS3, PLA2G16, AdPLA, H-REV107-1, HREV107, HREV107-3, MGC118754

描述:The HRAS like-suppressor 3 (HRASLS3) belongs to a class II tumor suppressor gene family and is involved in the regulation of differentiation and survival. This protein is expressed in several human tumors including ovarian carcinomas, lung carcinomas and may be involved in interferon-dependent cell death. Recombinant HRAS-like suppressor 3 was expressed in E.coli and was purified by conventional chromatography techniques.

配方:Liquid. In 20 mM Tris-HCl buffer (pH 8.0)

分子量:14.9 kDa (133 aa), confirmed by MALDI-TOF

序列:

MRAPIPEPKPGDLIEIFRPFYRHWAIIYVGDGYVVHLAPPSEVAGAGAASVMSALTDKAIIVKKELLYDVAGSDKYQVNNKH
DDKYSPLPCSKIIQRAEELVGQEVLYKLTSENCEHFVNELRYGVARSDQVRDV

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.