

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

HSBP1L1, 1-74aa, Human, His tag, E.coli

产品货号: TP02506

第三版

**别名:**Heat shock factor binding protein 1-like, FLJ10967, Heat shock factor-binding protein 1-like protein 1, MGC189743

**描述:**Heat shock factor binding protein 1-like, also known as HSBP1L1, belongs to a family of eukaryotic proteins identified as nucleotide exchange factors for HSP 70, which exhibit varying degrees of compartment and species specificity. This gene expression is regulated primarily at the transcription level. Recombinant human HSBP1L1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

**配方:**Liquid. In Phosphate buffer saline(pH 7.4) containing 10% glycerol

**分子量:**10.8 kDa (97aa) confirmed by MALDI-TOF

**序列:**

MGSSHHHHHHSSGLVPRGSHMGSM D VRGPEAPGGRALRDA AENLFQELQEHFQALTATLNLRMEEMGNRIEDLQKNV  
NDLMVQAGIENSIKEQMLKT

**纯度:**> 95% by HPLC

**浓度:**1 mg/ml (determined by Bradford assay)

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:**+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.