

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Interferon regulatory factor-3,1-112 aa, Human, Recombinant, E.coli

产品货号: TP02670

第三版

别名:IRF3

**描述:**Members of the Interferon regulatory factor (IRF) family regulate gene expression critical to immune response, hemopoiesis, and proliferation. IRF-3 is a member of the IRF family, and is distinct from other family members. Its transcriptional activity is regulated solely by posttranslational modifications. It plays a crucial role in activation of innate immunity and inflammation in response to viral infection.

**配方:**Liquid. In Phosphate-buffered saline (pH 7.4)

**分子量:**13.0kDa (112aa), confirmed by MALDI-TOF

**序列:**

**纯度:**> 95% by HPLC

**浓度:**1 mg/ml (determined by Bradford assay)

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.