

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Interferon- 1, Human, Recombinant, E.coli

产品货号: TP02672

第三版

别名:

描述: Interferon- 1 is one of the cytokines that are widely known to induce potent anti-viral activity . Interferon - 1 exerts a variety of other biological effects, including antitumor and immunomodulatory activities and are increasingly used clinically to treat a range of malignancies, myelodysplasias and autoimmune diseases. Recombinant human Interferon- 1 was expressed in E.coli and purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer.

配方: Liquid. In phosphate-buffered saline (PBS), pH 7.4

分子量: 19.5 kDa (167 amino acids)

序列:

纯度: > 95% by HPLC

浓度: 1mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.