

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Interleukin-6 receptor subunit beta, 23-617aa, Mouse, His tag, Insect cell

产品货号: TP02680

第三版

别名: IL6st, 5133400A03Rik, AA389424, BB405851, CD130, D13Ertd699e, gp130

**描述:** IL6st, also known as interleukin-6 receptor subunit beta, is a transmembrane protein which is the founding member of the class of all cytokine receptors. It forms one subunit of the type I cytokine receptor within the IL-6 receptor family. It is often referred to as the common gp130 subunit, and is important for signal transduction following cytokine engagement. As with other type I cytokine receptors, gp130 possesses a WSXWS amino acid motif that ensures correct protein folding and ligand binding. It interacts with Janus kinases to elicit an intracellular signal following receptor interaction with its ligand. Structurally, gp130 is composed of five fibronectin type-III domains and one immunoglobulin-like C2-type (immunoglobulin-like) domain in its extracellular portion. Recombinant mouse IL6st, fused to His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

**配方:** Liquid. In Phosphate Buffered Saline (pH 7.4) containing 10% glycerol.

**分子量:** 67.7kDa (603aa) 70-100kDa (SDS-PAGE under reducing conditions)

**序列:**

**纯度:** > 95% by HPLC

**浓度:** 1mg/ml (determined by Absorbance at 280nm)

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.