

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

JAZF1, 1-243aa, Human His tag, E.coli

产品货号: TP02706

第三版

别名:juxtaposed with another zinc finger protein 1, TIP27, ZNF802

描述:JAZF1 is a nuclear protein with three C2H2-type zinc fingers, and functions as a transcriptional repressor. Chromosomal aberrations involving this gene are associated with endometrial stromal tumors. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. Recombinant human JAZF1 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% glycerol

分子量:29.5kDa (266aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMGTGIAAASFFSNTCRFGGCGLHFPTLADLIEHIEDNHIDTDPRVLEKQELQQPTYVA
LSYINRFMTDAARREQESLKKKIQPKLSLTLSSSVSRGNVSTPPRHSSGSLTPPVTPPITPSSSFRSSTPTGSEYDEEEVDYEEESDS
DESWTTESAISSEAILSSMCMNGGEEKPFACVPVGCKKRYKNVNGIKYHAKNGHRTQIRVRKPFKCRGKSYKTAQGLRH
HTINFHPPVSAEIIIRKMQQ

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.