

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KISS1, 1-120aa, Human, Recombinant, E.coli

产品货号: TP02732

第三版

别名:KiSS-1 metastasis-suppressor, METASTIN, MGC39258

描述:KiSS-1 metastasis-suppressor (KISS1) is a metastasis suppressor gene that suppresses metastases of melanomas and breast carcinomas without affecting tumorigenicity. The metastasis suppressor properties may be mediated in part by cell cycle arrest and induction of apoptosis in malignant cells. Expression of KiSS1 has been localized to specific regions of the hypothalamus in many species and is regulated by gonadal steroids and across the estrous cycle.

配方:Liquid in 20mM Tris pH 8.0, 1mM DTT, 10% glycerol

分子量:12.6 kDa (120 aa) (Real molecular weight on SDS-PAGE will be shift up)

序列:

纯度:> 95% by HPLC

浓度:1 mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.