

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KRS (Lysyl-tRNA synthetase)

产品货号: TP02763

第三版

别名:

描述:Lysyl-tRNA synthetase (KRS) is a component of multi-ARS complex and catalyzes the ligation of lysine with lysyl-tRNA. KRS is secreted from cancer cell by TNF $\alpha$  and activates monocyte / macrophage. KRS is also a component of capsid proteins of HIV.

配方:1X PBS(pH7.2) 200mM NaCl 20% Glycerol

分子量:68kDa

序列:

纯度:> 95% by HPLC

浓度:0.50mg/ml

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.