

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LSM5, 1-91aa, Human, His tag, E.coli

产品货号: TP02852

第三版

别名: U6 snRNA-associated Sm-like protein LSm5 isoform a, LSM5 homolog, U6 small nuclear RNA associated (*S.cerevisiae*), YER146W

描述: Sm-like proteins were identified in a variety of organisms based on sequence homology with the Sm protein family. Sm-like proteins contain the Sm sequence motif, which consists of 2 regions separated by a linker of variable length that folds as a loop. The Sm-like proteins are thought to form a stable heteromer present in tri-snRNP particles, which are important for pre-mRNA splicing. Recombinant human LSM5 protein, fused to His-tag at N-terminus, was expressed in *E.coli* and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl(pH 8.5) containing 10% glycerol, 0.15M NaCl, 1mM DTT

分子量: 12.3 kDa (114aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSMANATTNPSQLLPLELVDKCIGSRIHIVMKSDKEIVGTLGFDFFVNMVLEDVTE
FEITPEGRRITKLDQILLNGNNITMLVPGGEGPEV

纯度: > 95% by HPLC

浓度: 1mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.