

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LY9, 48-454aa, Human, His-tag, Baculovirus

产品货号: TP02863

第三版

别名:T-lymphocyte surface antigen Ly-9 isoform , LY9, CD229, hly9, mLY9, SLAMF3, CD229Cell surface molecule Ly-9, Cell surface molecule Ly9hly9, LY 9, Ly9, LY9_HUMAN.

描述:LY9, as known as T-lymphocyte surface antigen Ly-9 isoform, is a type 1 transmembrane glycoprotein in the SLAM subgroup of the CD2 family. This protein consists of an extracellular domain with two Ig-like V-set and two Ig-like truncated C2-set domains. It is expressed on T and B cells, thymocytes, and more weakly on NK cells. Homophilic binding between LY9 molecules is mediated by the N-terminal Ig-like domain. Also, this protein knockout mice show minimally impaired immune responses, suggesting functional redundancy with other molecules. Recombinant human LY9, fused to His-tag at C-terminus, was expressed in insect cell and purified by using conventional chromatography techniques.

配方:Liquid. In Phosphate Buffered Saline (pH 7.4) containing 10% glycerol.

分子量:45.9kDa (416aa), 40-57kDa (SDS-PAGE under reducing conditions.)

序列:

ADLKDSAPTVVSGILGGSVTLPLNISVDTEIENVIWIGPKNALAFARPKENVTIMVKSYLGRDLITKWSYSLCISNLTNDAG
SYKAQINQRNFEVTTEEEFTLFVYEQLEPQVTMKSVKVSFSCNITLMCSVKGAEKSVLYSWTPREPHASESNGGSILTV
SRTPCDPLPYICTAQNPVSRSSLPVHVGQFCTDPGASRGGTTGETVVGVLGEPVTLPLALPACRDTEKVVWLFNTSIISK
EREEAATADPLIKSRDPYKNRVWVSSQDCSLKISQLKIEDAGPYHAYVCSEASSVTSMTHVTLIIYRRLRKPKITWLSLRHSED
GICRISLTCSVEDGGNTVMYTWTPQLQKEAVVSQGESHLNVSWRSENHPNLTCTASNPVSRSSHQFLSENICSGPERNTKH
HHHHH

纯度:> 95% by HPLC

浓度:0.25mg/ml (determined by Absorbance at 280nm)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储:+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.