

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

LYVE1, 25-235aa, Human, His tagged, Recombinant, Insect cell

产品货号: TP02870

---

第三版

**别名:**Lymphatic vessel endothelial hyaluronic acid receptor 1, CRSBP-1, HAR

**描述:**Lymphatic vascular endothelial hyaluronan receptor-1 (LYVE-1), a selective marker of the lymphatic endothelium, is a surface endocytic receptor for hyaluronan (HA), which is an extracellular glycosaminoglycan involved in cell adhesion and migration. LYVE1 is also known as XLKD1 and may function in lymphatic hyaluronan transport and have a role in tumor metastasis.

**配方:**Liquid in 20mM Tris pH7.5, 10% Glycerol

**分子量:**24.8 kDa (229 aa)

**序列:**

**纯度:**> 95% by HPLC

**浓度:**1 mg/ml

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.