



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MAK16, 1-300aa, Human, His tag, E.coli

产品货号: TP02889

第三版

别名: MAK16 homolog, MAK16L, RBM13

描述: MAK16 belongs to the MAK16 family. Recombinant human MAK16 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In Phosphate buffered saline (pH7.4), 30% glycerol, 1mM DTT

分子量: 37.8 kDa (323aa) confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMQSDDVWDTLGNKQFCSFKIRTKTQSFCRNEYSLTGLCNRSSCPLANSQYATIKEE
KGQCYLYMKVIERAAFPRRLWERVRLSKNYEKALEQIDENLIYWPRFIRHKCKQRFTKITQYLIRIRKLT LKRQRKLVPLSKK
VERREKRREEKALIAAQLDNAIEKELLERLQDQTYGDIYNFPIHAFDKALEQQEAESDSSDTEEKDDDDDDDEEDVGKREFV
EDGEVDES DISDFEDMDKLDASDDEDQDGKSSSEEEEEKALSAKHKGKMP LRGPLQRKRAYVEIEYEQETEPVAKAKTT

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.