

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MEIS3, 1-358aa Human, His tag, E.coli

产品货号: TP02944

第三版

别名: Homeobox protein Meis 3, MRG2, DKFZp547H236, Homeobox protein Meis3, MRG2, Meis1-related protein 2, Meis homeobox 3

描述: MEIS3 belongs to the TALE/MEIS homeobox family and contains 1 homeobox DNA-binding domain. Recombinant human MEIS3 protein, fused to His-tag at N-terminus, was expressed in E.coli .

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

分子量: 41.4 kDa (381 aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMARRYDELPHYPGIVDGPAAALASFPETVPAVPGPYGPHRPPQPLPPGLDSDGLKRE
KDEIYGHPLFPLLALVFEEKCELATCSPRDGAGAGLGTTPGGDVCSSDSFNEDIAAFAKQVRSERPLFSSNPELDNLMIQAIQV
LRFHLLLEKKGKMPIDLVIEDRDGGCREDFEDYPASCPSLPDQNNMWIRDHEDSGSVHLGTPGPSSGGLASQSGDNSSDQ
GDGLDTSVASPSSGGEDEDLDQERRRNKKGIFPKVATNIMRAWLWFQHLSPYPSEEQKQLAQDTGLTILQVNNWFINA
RRRIVQPMIDQSNRTGQGAAFSPEGQPIGGYTETQPHVAVRPPGSVGM SLNLEGEWHYL

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.