

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

METTL21A, 93-218aa, Human, His tag, E.coli

产品货号: TP02957

第三版

别名: Methyltransferase-like protein 21A, FAM119A, HCA557b

描述: METTL21A belongs to the methyltransferase superfamily and is probable methyltransferase. Recombinant human METTL21A protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% glycerol.

分子量: 17 kDa (149aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSTDRKVALEFLKSNVQANLPPHIQTKTVVKELTWGQNLGSFSPGEFDLILGADIILE
ETFTDLLQLEHLCSNHSVILLACRIRYERDNNFLAMLERQFTVRKVHYDPEKDVHIYEAQKRQKEDL

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.