



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MFAP3, 19-147aa, Human, His tag, E.coli

产品货号: TP02959

第三版

**别名:** Microfibril-associated glycoprotein 3 isoform 1, microfibrillar-associated protein 3, Microfibril associated glycoprotein 3

**描述:** MFAP3 contains 1 Ig-like C2-type (immunoglobulin-like) domain. This protein is a component of the elastin-associated microfibrils. Recombinant human MFAP3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

**配方:** Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT.

**分子量:** 16.7 kDa (152aa) confirmed by MALDI-TOF

**序列:**

MGSSHHHHHHSSGLVPRGSHMGSAFVLEDVDFDQMVSLANRSSYNASFPSSFELSASSHSDDDVIHAKGTSVSIECLLTA  
SHYEDVHWHNSKGGQLDGRSRGGKWLVSDFLNITNVAFDDRGLYTCFVTSPIRASYSVTLRVIFTSGDM

**纯度:** > 95% by HPLC

**浓度:** 0.5 mg/ml (determined by Bradford assay)

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.