



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MIA, 25-131aa, Human, His tag, E.coli

产品货号: TP02969

第三版

别名: Melanoma inhibitory activity, CD-RAP

描述: MIA belongs to the MIA/OTOR family and contains 1 SH3 domain. This protein elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas. Recombinant human MIA protein, fused to His-tag at N-terminus, was expressed in E.coli

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M Urea

分子量: 14.4 kDa (128aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGPMPLADRKLCADQECSPISMAVALQDYMAPDCRFLTIHRGQVVVYVFSKLGKRG  
RLFWGGSVQGDYYGDLAARLGYFPSSIVREDQTLKPGKVDVKTDKWDYFCQ

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.