

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MOSPD1, 1-158aa, Human, His tag, E.coli

产品货号: TP03009

第三版

别名: Motile sperm domain containing 1, DJ473B4

描述: MOSPD1 contains 1 MSP domain. Recombinant human MOSPD1 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方: Liquid. In Phosphate buffered saline (pH7.4) containing 0.15M NaCl, 50% glycerol, 1mM DTT

分子量: 20.6kDa (181aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMHQKQRPELVEGNLPVVFVPTTELIFYADDQSTHKQVLTLYNPYEFALKFKVLCTT
PNKYVVVDAAGAVKPCQCCVDIVIRHRDVRSCYGVIDKFRLQVSEQSQRKALGRKEVVATLLPSAKEQQKEEEEKRLKE
HLTESLFFEQSFQPENRAVSSGP

纯度: > 95% by HPLC

浓度: 0.5mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.