

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MRFAP1L1, 1-127aa Human, His tag, E.coli

产品货号: TP03017

第三版

别名:MORF4 family-associated protein 1-like, Morf4 family associated protein 1-like 1, PP784

描述:MRFAP1L1 belongs to the MORF4 family-associated protein family. Recombinant human MRFAP1L1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

分子量:17.2 kDa (150aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMRPLDIDEVEAPEEVEVLEPEEDFEQFLLPVINEMREDIASLIREHGRAYLRTRSKLW
EMDNMLIQIKTQVEASEESALNHVQHPSGEADERVSELCEKAEKAKEIAKMAEMLVELVWRIERSESS

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.