

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MRPL1, 51-325aa, Human, His tag, E.coli

产品货号: TP03021

第三版

别名:39S ribosomal protein L1, mitochondrial precursor , BM022, L1MT, MRP-L1

描述:39S ribosomal protein L1, mitochondrial precursor, also known as MRPL1, is a 338 amino acid mitochondrial protein that exists as a component of the 39S ribosomal subunit and works in conjunction with other MRPs to mediate protein synthesis. Mitochondrial ribosomes consist of a large 39S subunit and a small 28S subunit, both of which are comprised of multiple mitochondrial ribosomal proteins (MRPs) that are encoded by nuclear genes and are essential for protein synthesis within mitochondria. MRPL1 exists as two isoforms produced by alternative splicing. Isoform one of MRPL1 is ubiquitously expressed while isoform two is specifically expressed in heart. Recombinant human MRPL1 protein to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 30% glycerol, 1mM DTT

分子量:33.8kDa (298aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHSSGLVPRGSHMGSKKTKKGAKEKTPDEKKDEIEKIKAYPYMEGEPEDDVYLKRLYPRQIYEVEKAVHLL
KKFQILDFTSPKQSVYLDLTLDMALGKKKNVEPFTSVLSLPYPFASEINKVAVFTENASEVKIAEENGAAFAGGTSLIQKIWD
DEIVADFYVAVPEIMPELNRLRKKLNKKYPKLSRNSIGRDIPKMLELFKNGHEIKVDEERENFLQTKIATLDMSSDQIAANL
QAVINEVCRHRPLNLGPFVVRAFLRSSTSEGLLLKIDPLLPKEVKNEESEKEDA

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.