

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MRPS25, 1-173aa, Human, His tag, E.coli

产品货号: TP03028

第三版

别名:28S ribosomal protein S25, mitochondrial, mitochondrial ribosomal protein S25, MRP-S25, RPMS25

描述:MRPS25 is a 28S subunit protein that belongs to the ribosomal protein S5P family. A pseudogene corresponding to this gene is found on chromosome 4. Recombinant human MRPS25 protein, fused to His-tag at N-terminus, was expressed in E.coli .

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing, 10% glycerol, 0.4M Urea

分子量:22.5kDa (196aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMMPKGRFPIRRTLQYLSQGNVVKDSVKVMTVNYNTHGELGEGARKFVFFNIPQI
QYKNPWVQIMMFKNMTPSPFLRFYLDSEQLVDVETKSNKEIMEHIRKILGKNEETLREEEEEKKQLSHPANFGPRKYCL
RECICEVEGQVPCPSLVPLPKEMRGKYKAALKADAQD

纯度:> 95% by HPLC

浓度:1 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.