

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MUCL1, 21-90aa Human, His tag, E.coli

产品货号: TP03046

第三版

别名: Mucin-like protein 1, SBEM

描述: MUCL1 may play a role as marker for the diagnosis of metastatic breast cancer. This protein is expressed in mammary, salivary glands and prostate. Recombinant human MUCL1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In Phosphate buffer saline (pH 7.4) containing 10% glycerol

分子量: 9.3 kDa (93aa) confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSNPTTAAPADTYPATGPPADDEAPDAETTAATTATTAAPTTATTAASTTARKDIPVL
PKWVGDLPNGRVCP

纯度: > 95% by HPLC

浓度: 0.25 mg/ml (determined by BCA assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.