

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

NDRG3, 1-363aa, Human, His tag, E.coli

产品货号: TP03121

第三版

别名:Protein NDRG3, NDRG family member 3, N-myc downstream-regulated gene 3 protein

描述:NDRG3 belongs to the NDRG family and is highly expressed in brain. Recombinant human NDRG3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques

配方:Liquid. In Phosphate buffered saline (pH7.4) containing 20% glycerol, 1mM DTT

分子量:42.4 kDa (386aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMDELQDVQLTEIKPLLNDKEHDIETTHGVVHVTIRGLPKGNRPVILTYHDIGLNH
KSCFNAAFFNFEDMQEITQHFAVCHVDAPGQQEGAPSFPTGYQYPTMDELAEMLPVLTHLSLKSIIIGIGVGAGAYILSRFAL
NHPELVEGLVLINVDPCAAGWIDWAASKLSGLTTNVVDIILAHHFGQEELQANLDLIQTYRMHIAQDINQDNLQFLNSY
NGRRDLEIERPILGQNDNKSRTLKSTLLVVGDNSPA VEAVVECNRLNPINTLLKMADCGGLPQVVQPGKLTEAFKYFL
QGMGYIPASMTRLARSRTHSTSSSLGSGESPFSSVTSNQSDGTQESCESPDVDRHQTM EVSC

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.