

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

OCIAD2, 1-154aa, Human, His tag, E.coli

产品货号: TP03243

第三版

别名: OCIA domain-containing protein 2 isoform 1, MGC45416, OCIA domain-containing protein 2, Ovarian carcinoma immunoreactive antigen-like protein

描述: OCIAD2 is a protein that localizes to endosomes and contains one OCIA (Ovarian carcinoma immunoreactive antigen) domain. Recombinant human OCIAD2 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% glycerol

分子量: 19.3kDa (177aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMASASARGNQDKDAHFPSPKQSLLFCKPKSLHIHRAEISKIMRECQEESFWKRALP
FSLVSMMLVTQGLVYQGYLAANSRFGSLPKVALAGLLGFLGKVSIGVCQSKFHFFEDQLRGAGFGPQHNRHCLLTCEEC
KIKHGLSEKGDSPSAS

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.