

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

OXLD1, 46-147aa, Human, His tag, E.coli

产品货号: TP03273

第三版

别名:Oxidoreductase-like domain-containing protein 1, C17orf90

描述:OXLD1 contains 1 Oxidoreductase-like domain. Recombinant human OXLD1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 50% glycerol, 1mM DTT

分子量:13.4 kDa (125aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSPGAQAPDGRRKFGTDHVEVGSQAGADGTRPPKASLPPELQPPTNCCMSGCPNCV
WVEYADRLLQHFQDGGERALAALEEHVADENLKAFLRMEIRLHTRCGG

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.