

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Parathyroid Hormone Recombinant, E.coli

产品货号: TP03290

第三版

别名: PTH, PTH1.

描述: Human parathyroid hormone (hPTH) is an 84 amino acid residue peptide and one of the main regulators in maintenance of calcium homeostasis. It acts primarily on kidney and bone cells stimulating calcium back resorption or calcium mobilization, respectively. The whole mechanism is still under discussion, but low-dose hPTH triggers cyclic AMP-dependent protein kinase in some populations of bone cells bearing PTH receptors, which stimulates the proliferation of osteoblasts.

配方: Liquid. In phosphate-buffered Saline (pH 7.4) containing 2 mM EDTA.

分子量: 9.55 kDa (85 aa), confirmed by MALDI-TOF

序列:

MSVSEIQLMHNLGKHLNSMERVEWLRKKLQDVHNFVALGAPLAPRDAGSQRPRKKEDNVLVESHEKSLGEADKADVNVLTAKAKSQ

纯度: > 95% by HPLC

浓度: 1 mg/ml (determined by Bradford assay)

内毒素: < 1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.