

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PCP4L1, 1-68aa, Human, His tag, E.coli

产品货号: TP03311

第三版

别名: Purkinje cell protein 4-like protein 1, IQM1

描述: PCP4L1, also known as Purkinje cell protein 4-like protein 1, belongs to the PCP4 family and contains 1 IQ domain. Recombinant human PCP4L1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol

分子量: 9.9kDa (91aa), confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMSELNKTSPATNQAAGQEEKGKAGNVKKAEEEEIDIDLTAPETEKAALAIQGKF  
RRFQKRKKDPSS

纯度: > 95% by HPLC

浓度: 0.5 mg/ml (determined by BCA assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.