



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PDGFD, 250-370aa, Human, His tag, E.coli

产品货号: TP03324

第三版

别名: Platelet-derived growth factor D, IEGF, MSTP036, SCDGF-B, SCDGFB.

描述: Platelet-derived growth factor D, also known as PDGFD, plays an important role in wound healing. This protein induces macrophage recruitment, increased interstitial pressure, and blood vessel maturation during angiogenesis. This protein can initiate events that lead to a mesangial proliferative glomerulonephritis, including influx of monocytes and macrophages and production of extracellular matrix.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

分子量: 16.6 kDa (146aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSHMSYHDRKSKVDLRLNDDAKRYSCTPRNYSVNIREELKLANVVFPRCLLVQRC
GGNCGCGTVNWRSCTCNSGKTVKKYHEVLQFEPGHIKRRGRAKTMALVDIQLDHERCDCICSSRPPR

纯度: > 95% by HPLC

浓度: 1 mg/ml

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.