

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PIH1D2, 1-315aa Human, His tag, E.coli

产品货号: TP03406

第三版

别名:PIH1 domain-containing protein 2, PIH1 domain containing 2

描述:PIH1D2 belongs to the PIH1 family. Recombinant human PIH1D2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid, In Phosphate buffered saline (pH7.4) containing 20% glycerol, 1mM DTT

分子量:38.3 kDa (338aa) confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMETSSKGLLTQVTQFWNLLDDLAQSDPEGYEFIQQLKEGKQLCAAPEPQLCLQ  
TRILKPKEKILFINLCQWTRIPAPQSTTHPVPLTVGKPEDTTEISDAYTVIDVAYNPDVLAHAAEKDQVKKNQLIQMAMKCI  
EEKFQFTLSHSYHITKFRIKGSIQRMKQNLMGITDSDILREKMRRELTLGQIRSSTMSNPDPHFPQLLLPKDQVSGKAVCLIE  
EISSTEIQVEMKMPAYELKIVHDHSEKPLKIELKVELPGINSVSLCDLSVSEDDLLIEVSEKYRLHLNLPKLIDEMTTAKFIKE  
KSTLIITMPLV

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.