

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PLAC8, 1-115aa, Human, His tag, E.coli

产品货号: TP03425

第三版

别名: Placenta-specific 8, C15, onzin

描述: PLAC8 belongs to the cornifelin family. This protein is expressed at high levels in plasmacytoid dendritic cells, spleen, lymph nodes, peripheral blood leukocytes, and bone marrow. Recombinant human PLAC8 protein, fused to His-tag at N-terminus, was expressed in E.coli.

配方: Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M urea

分子量: 14.9 kDa (138aa)

序列:

MGSSHHHHHHSSGLVPRGSHMGSMQAQAPVVVVVTQPGVGPAPQNSNWQTGMCDGFSDCGVCLCGTFCFPCLGC
QVAADMNECCLCGTSVAMRTLYRTRYGIPGSICDDYMATLCCPHCTLCQIKRDINRRRAMRTF

纯度: > 95% by HPLC

浓度: 1mg/ml (determined by Bradford assay)

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.