



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PLCXD3, 1-321aa, Human, His tag, E.coli

产品货号: TP03434

第三版

别名:PI-PLC X domain-containing protein 3, phosphatidylinositol-specific phospholipase C X domain containing 3, phosphatidylinositol-specific phospholipase C, X domain containing 3

描述:PLCXD3 contains 1 PI-PLC X-box domain. Recombinant human PLCXD3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

分子量:38.7 kDa (344aa) confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMASQKNEKLADWMATLPESMHSIPLTNLAIPGSHDSFSFYIDEASPVGPEQPET
VQNFVSVFGTVAKKLMRKWLATQTMNFTGQLGAGIRYFDLRISTKPRDPDNELYFAHGLFSAKVNEGLEEINAFITDHH
KEVVFLDFNHFYGMQKYHHEKLVQMLKDIYGNKMCPAIFAQEVSLKYLWEKDYQVLVIFYHSPVALEVPFLWPGQMMP
APWANTTDPEKLIQFLQASITERKKGSFFISQVVLTPKASTVVKGVASGLRETITERALPAMMQWVRTQKPGESGINIVTA
DFVELGDFISTVIKLNYYVFDEGEANT

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储:+4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.