

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PLEKHF2, 1-249aa Human, His tag, E.coli

产品货号: TP03436

第三版

别名:pleckstrin homology domain-containing family F member, EAPF, PHAFIN2, ZFYVE18

描述:PLEKHF2 contains 1 FYVE-type zinc finger and 1 PH domain. Recombinant human PLEKHF2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

分子量:30.2kDa (272aa), confirmed by MALDI-TOF

序列:

MGSSHHHHHHSSGLVPRGSHMGSMVDRLANSEANTRRISIVENCFGAAGQPLTIPGRVLI GEGVLT KLCRKKPKARQFFL
FNDILVYGNIVIQKKKYNKQHIIPLNVTIDSIKDEGDLRNGWLIKTPTKSAVYAATATEKSEWMNHINKCVTDLLSKSG
KTPSNEHA AVWVPDSEATVCMRCQKAKFTPVNRRHHCRKCGFVVCGPCSEKRFLLP SQSSKPV RICDFCYDLLSAGDMA
TCQPARDSDSYSQSLKSPLNDMSDDDDDDDDSSD

纯度:> 95% by HPLC

浓度:0.25 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.