

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

POLR3H, 1-204aa

产品货号: TP03469

第三版

别名:DNA-directed RNA polymerase III subunit RPC8, RPC22.9, RPC8.

描述:POLR3H, also as known as A1665 and RPC8, belongs to the eukaryotic RPB7/RPC8 RNA polymerase subunit family. POLR3H is DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. It plays a key role in sensing and limiting infection by intracellular bacteria and DNA viruses. POLR3H acts as nuclear and cytosolic DNA sensor involved in innate immune response. Recombinant human POLR3H protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

配方:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT

分子量:25.3kDa (227aa) confirmed by MALDI-TOF (Molecular size on SDS-PAGE will appear higher)

序列:

MGSSHHHHHSSGLVPRGSHMGSMFVLVEMVDTVRIPPWQFERKLNDISIAEELNKKLANKVVYINVGLCICLFDITKLED
AYVFPDGDASHTKVFHRCVVFHFPFLDEILIGKIKGCSPEGVHVS LGFFDDILIPPELQQPAKFDEAEQVWVWEYETEEGA
HDLYMDTGEEIRFRVVDSEFVDTSPGSSADATTSSEELPKKEAPYTLVGSISEPGLGLLSWWTSN

纯度:> 95% by HPLC

浓度:0.5 mg/ml (determined by Bradford assay)

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储:+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.