



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Protein Marker-1

产品货号: TP03541

第三版

别名:n/a

描述:Mid-Range Protein Marker-1 consists of Seven precisely sized recombinant proteins (M.W. 70, 57, 40, 28, 18, 13.5, 8.5KDa). The concentration of each protein is 0.1ug/ul except for 8.5KDa (0.2ug/ul). Usage Note: Use 5ul per lane and you need not heat before use.

配方:12mM Tris-HCl pH6.8, 5% glycerol, 0.4% SDS 2.88 mM 2-mercaptoethanol, 0.02% bromophenol blue

分子量:34

序列:

纯度:> 95% by HPLC

浓度:n/a

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.