

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Protein Marker-2, Broad range

产品货号: TP03542

第三版

别名:

描述: ATGen Broad-Range Protein Marker-2 consists of eight precisely sized recombinant proteins (M.W. 150, 100, 70, 50, 35, 25, 20, 15 KDa). The concentration of each protein is about  $0.1 \mu\text{g}/\mu\text{l}$ . For research use only, not for use in diagnostic procedures.

配方: Please see datasheet.

分子量:

序列:

纯度:  $> 95\%$  by HPLC

浓度: Please see datasheet.

内毒素:  $< 1.0$  EU per 1  $\mu\text{g}$  of protein (determined by LAL method)

存储:  $+4^\circ\text{C}$  保存 (1-2 周). 长期保存在  $-20^\circ\text{C}$  或者  $-70^\circ\text{C}$ . 避免反复冻融.